

# IGFBP7 ELISA Kit (Human) (OKBB00381) Lot# KE0897

# **Instruction for Use**

For the quantitative measurement of Human IGFBP7 in cell culture supernatants, serum, plasma (heparin, EDTA) and urine.

Variation between lots can occur. Refer to the manual provided with the kit.

This product is intended for research use only.



## **Table of Contents**

1.	Background2
2.	Assay Summary3
3.	Storage and Stability3
4.	Kit Components3
5.	Precautions4
6.	Required Materials Not Supplied4
7.	Technical Application Tips4
8.	Reagent Preparation5
9.	Sample Preparation7
10.	Assay Procedure8
11.	Calculation of Results9
12.	Typical Expected Data9
13.	Technical Resources



## 1. Background

### **Principle**

Aviva Systems Biology IGFBP7 ELISA Kit (Human) (OKBB00381) is based on standard sandwich enzymelinked immuno-sorbent assay technology. A mouse monoclonal antibody specific for IGFBP7 has been precoated onto 96-wellplate (12 x 8 Well Strips). Standards (NSO; D30-L282) and test samples are added to the wells and incubated. After washing, a biotinylated polyclonal goat detector antibody specific for IGFBP7 is added, incubated and followed by washing. Avidin-Biotin-Peroxidase Complex is then added, incubated and unbound conjugate is washed away. An enzymatic reaction is visualized through the addition of TMB substrate which is catalyzed by HRP to produce a blue color product that changes yellow after adding acidic stop solution. The density of yellow coloration read by absorbance at 450 nm and is quantitatively proportional to the amount of sample Human IGFBP7 captured in well.

## **Background**


Insulin-like growth factor-binding protein 7 is a protein that in humans is encoded by the IGFBP7 gene. The gene was mapped to chromosome 4q12. The major function of the protein is the regulation of availability of insulin-like growth factors (IGFs) in tissue as well as in modulating IGF binding to its receptors. IGFBP7 binds to IGF with high affinity. It also stimulates cell adhesion. The protein is implicated in some cancers.

### **General Specifications**

General Specificiations						
Range 0.625 ng/mL -40 ng/mL						
Sensitivity < 0.020 ng/mL (Derived by linear regression of OD <sub>450</sub> of the Mean Blank + 2xSD)						
	Natural and recombinant Human IGFBP7					
	UniProt ID: Q16270					
Specificity	Gene ID: 602867					
	<u>Target Alias</u> : AGM, FSTL2, IBP-7, IGF-binding protein 7, IGFBP-7, IGFBP-7v, IGFBPRP1, IGFBP-rP1, Insulin-like growth factor-binding protein 7, MAC25, MAC25 protein, PGI2-stimulating factor, Prostacyclin-stimulating factor, PSF, TAF, Tumor-derived adhesion factor					
Cross-Reactivity No detectable cross-reactivity with other relevant proteins						



## 2. Assay Summary



## 3. Storage and Stability

• Upon receipt store kit at -20°C for 6 months. Avoid multiple freeze/thaw cycles.

## 4. Kit Components

• The following reagents are the provided contents of the kit.

Description	Quantity	Storage Conditions	
96-Well plate Pre-coated with Anti-Human IGFBP7 Antibody	1 (12 x 8 Well Strip)		
Lyophilized Recombinant Human IGFBP7 standard	2 x 40 ng		
100X Biotinylated Anti-Human IGFBP7 Antibody	1 x 130 μL		
100X Avidin-Biotin-Peroxidase Complex (ABC)	1 x 130 μL		
Sample Diluent Buffer	1 x 30 mL	Store at -20°C for 6 months	
Antibody Diluent Buffer	1 x 12 mL		
ABC Diluent Buffer	1 x 12 mL		
TMB Color Developing Agent	1 x 10 mL		
TMB Stop Solution	1 x 10 mL		
10X Wash Buffer	1 x 30 mL		



#### 5. Precautions

- Read instructions fully prior to beginning use of the assay kit.
- Any deviations or modifications from the described method or use of other reagents could result in a reduction of performance.
- Reduce exposure to potentially harmful substances by wearing personal protective lab equipment including lab coats, gloves and glasses.
- For information on hazardous substances included in the kit please refer to the Material Safety Data Sheet (MSDS).

## 6. Required Materials Not Supplied

- Microplate reader capable of reading absorbance at 450 nm.
- Automated plate washer (optional).
- Pipettes capable of precisely dispensing 0.5 µL through 1 mL volumes of aqueous solutions.
- Pipettes or volumetric glassware capable of precisely measuring 1 mL through 100 mL of aqueous solutions.
- New, clean tubes and/or micro-centrifuge tubes for the preparation of standards or samples.
- Absorbent paper or paper toweling.
- Distilled or deionized ultrapure water.

## 7. Technical Application Tips

- Do not mix or substitute components from other kits.
- To ensure the validity of experimental operation, it is recommended that pilot experiments using standards and a small selection of sample dilutions to ensure optimal dilution range for quantitation.
- Samples exhibiting OD measurements higher than the highest standard should be diluted further in the appropriate sample dilution buffers.
- Inspect all reagents prior to use. Components should contain no particulates or cloudiness and should be colorless.
- Prior to using the kit, briefly spin component tubes to collect all reagent at the bottom.
- Replicate wells are recommended for standards and samples.
- Cover microplate while incubating to prevent evaporation.
- •Do not allow the microplate wells dry at any point during the assay procedure.
- Do not reuse tips or tube to prevent cross contamination.
- •Avoid causing bubbles or foaming when pipetting, mixing or reconstituting.
- Completely remove of all liquids when washing to prevent cross contamination.
- Prepare reagents immediately prior to use and do not store, with the exception of the top standard.
- Equilibrate all materials to ambient room temperature prior to use (standards exception).
- For optimal results in inter- intra-assay consistency, equilibrate all materials to room temperature prior to performing assay (standards exception) and perform all incubations at 37°C.
- Pipetting less than 1 µL is not recommended for optimal assay accuracy.
- Once the procedure has been started, all steps should be completed without interruption. Ensure that all reagents, materials and devices are ready at the appropriate time.
- Incubation times will affect results. All wells should be handled in the same sequential order and time intervals for optimal results.
- Samples containing precipitates or fibrin strands or which are hemolytic of lipemic might cause inaccurate results due to interfering factors.




## 8. Reagent Preparation

Equilibrate all materials to room temperature prior to use and use immediately.

### 8.1 Human IGFBP7 Assay Standards

- 8.1.1 Prepare the IGFBP7 standards no greater than 2 hours prior to performing experiment. Standards should be held on ice until use in the experiment.
- 8.1.2 Reconstitute one of the provided 40 ng Lyophilized Recombinant Human IGFBP7 Standard for each experiment. Prepare a stock 40 ng/mL Human IGFBP7 Standard by reconstituting one tube of Lyophilized Recombinant Human IGFBP7 Standard as follows:
  - 8.1.2.1 Gently spin or tap the vial to collect all material at the bottom.
  - 8.1.2.2 Add 1 mL of **Sample Diluent Buffer** to the vial.
  - 8.1.2.3 Seal then mix gently and thoroughly.
  - 8.1.2.4 Leave the vial to sit at ambient temperature for 10 minutes.
- 8.1.3 Prepare a set of seven serially diluted standards as follows:
  - 8.1.3.1 Label tubes with numbers 2 8.
  - 8.1.3.2 Add 300  $\mu$ L of **Sample Diluent Buffer** to Tube #'s 2 8.
  - 8.1.3.3 Use the reconstituted 40 ng/mL Human IGFBP7 as **Standard #1**.
  - 8.1.3.4 Prepare **Standard #2** by adding 300 μL of **40 ng/mL Standard #1** from Tube #1 to Tube #2. Mix gently and thoroughly.
  - 8.1.3.5 Prepare **Standard #3** by adding 300  $\mu$ L of **Standard #2** from Tube #2 to Tube #3. Mix gently and thoroughly.
  - 8.1.3.6 Prepare further serial dilutions through Tube #7. Reference the table below as a guide for serial dilution scheme.
  - 8.1.3.7 Tube #8 is a blank standard (only **Sample Diluent Buffer**), which should be included with every experiment.

Standard Number (Tube)	Sample To Dilute	Volume standard (µL)	Volume Sample Diluent Buffer (μL)	Total Volume (μL)	Final Concentration
Tube #1	40 ng/mL of Human IGFBP7 standard	NA	NA	1,000	40 ng/mL
Tube #2	40 ng/mL	300	300	600	20 ng/mL
Tube #3	20 ng/mL	300	300	600	10 ng/mL
Tube #4	10 ng/mL	300	300	600	5 ng/mL
Tube #5	5 ng/mL	300	300	600	2.5 ng/mL
Tube #6	2.5 ng/mL	300	300	600	1.25 ng/mL
Tube #7	1.25 ng/mL	300	300	600	0.625 ng/mL
Tube #8	NA	0	300	300	0.0 (Blank)





#### 8.2 1X Biotinylated Anti-Human IGFBP7 Antibody

- 8.2.1 Prepare the **1X Biotinylated Anti-Human IGFBP7 Antibody** immediately prior to use by diluting the **100X Biotinylated Anti-Human IGFBP7 Antibody** 1:100 with **Antibody Diluent Buffer**.
- 8.2.2 For each well to be used in the experiment prepare 120 μL by adding 1.2 μL of **100X** Biotinylated Anti-Human IGFBP7 Antibody to 118.8 μL Antibody Diluent Buffer.
- 8.2.3 Mix thoroughly and gently. Hold no longer than 2 hours prior to using in procedure.

### 8.3 1X Avidin-Biotin-Peroxidase Complex (ABC)

- 8.3.1 Prepare the **1X Avidin-Biotin-Peroxidase Complex (ABC)** immediately prior to use by diluting the **100X Avidin-Biotin-Peroxidase Complex (ABC) 1:100** with **ABC Dilution Buffer**.
- 8.3.2 For each well to be used in the experiment prepare 120 μL, by adding 1.2 μL of **100X Avidin-Biotin-Peroxidase Complex (ABC)** to 118.8 μL **ABC Dilution Buffer**.
- 8.3.3 Mix thoroughly and gently. Hold no longer than 2 hours prior to using in procedure.

### 8.4 1X Wash Buffer

- 8.4.1 Add 270 mL of ultra-pure water to a clean > 500 mL bottle or other vessel.
- 8.4.2 Add the entire 30 mL contents of the **10X Wash Buffer** bottle to the water.
- 8.4.3 Seal and mix gently by inversion. Avoid foaming or bubbles.
- 8.4.4 Store the **1X Wash Buffer at** room temperature until ready to use in the procedure. Store the prepared **1X Wash Buffer** at 4°C for no longer than 1 week. Do not freeze.

#### 8.5 Microplate Preparation

- Micro-plates are provided ready to use and do not require rinsing or blocking.
- Unused well strips should be returned to the original packaging, sealed and stored at 4°C.
- Equilibrate microplates to ambient temperatures prior to opening to reduce potential condensation.



## 9. Sample Preparation

#### 9.1 Sample Preparation and Storage

- Store samples to be assayed at 2-8°C for 24 hours prior being assayed.
- For long term storage, aliquot and freeze samples at -20°C. Avoid repeated freeze-thaw cycles.
- Clear samples by centrifugation as follows:
  - **Cell culture supernatants** Remove particulates by centrifugation, assay immediately or aliquot and store samples at -20°C.
  - **Serum** Allow the serum to clot in a serum separator tube (about 4 hours) at room temperature. Centrifuge at approximately 1,000 x g for 15 min. Analyze the serum immediately or aliquot and store frozen at -70°C.
  - Plasma Collect plasma using heparin or EDTA as an anticoagulant. Centrifuge for 15 min at 1,000 x g within 30 min of collection. Analyze immediately or aliquot and store frozen at -70°C. Citrate is not recommended as the anticoagulant.
  - **Urine** -Aseptically collect the first urine of the day, micturate directly into a sterile container. Remove particular impurities by centrifugation, assay immediately or aliquot and store samples at -20°C.

#### 9.2 Sample Dilution

Target protein concentration must be estimated and appropriate sample dilution selected such that the final target protein concentration falls near the middle of the assay linear dynamic range.

- Prepare 150 µL sample for each replicate to be assayed.
- Dilute samples with Sample Diluent Buffer.
- · Mix diluted samples gently and thoroughly.
- Pipetting less than 2 µL is not recommended for optimal assay accuracy.
- Refer to the following table for recommended sample dilution guidelines based on the dynamic range of this kit:

Estimated Sample T	arget Concentration	Dilution Level	Sample Volume For Two Replicates	Sample Diluent Buffer For Two Replicates	
High Concentration	400-4,000 ng/mL	1:100	3 µL	297 μL	
Medium Concentration	40-400 ng/mL	1:10	25 µL	225 μL	
Low Concentration	625-40 ng/mL	1:2	100 μL	100 μL	
Very Low Concentration	≤ 0.625 ng/mL	1:2 or No Dilution	-	-	



## 10. Assay Procedure

- The ABC Working Solution, TMB Color Developing Agent and TMB Stop Solution must be kept warm at 37°C for 30 minutes prior to use. Equilibrate all other reagents and materials to ambient room temperature prior to use in the procedure.
- To control for small potential variations in micro well-plate and day to day ambient temperature fluctuations, equilibrate all reagents prior to use and perform all incubation steps at to 37°C for optimal consistency and reproducibility.
- **10.1** Add 100 μL of serially titrated standards, diluted samples or blank into wells of the pre-coated well plate. At least two replicates of each standard, sample or blank is recommended.
- **10.2** Cover the plate with the plate sealer and incubate at 37°C for 90 minutes.
- **10.3** Remove the plate sealer and discard the liquid in the wells by rigorously flicking into an acceptable waste receptacle.
- **10.4** Gently blot any remaining liquid from the wells by tapping inverted on the benchtop onto paper toweling. Do not allow the wells to completely dry at any time.
- 10.5 Add 100 µL of prepared 1X Biotinylated Anti-Human IGFBP7 Antibody to each well.
- **10.6** Cover with the plate sealer and incubate at 37°C for 60 minutes.
- 10.7 Wash plate 3 times with 1X Wash Buffer as follows:
  - 10.7.1 Discard the liquid in the wells by rigorously flicking into an acceptable waste receptacle.
  - 10.7.2 Gently blot any remaining liquid from the wells by tapping inverted on the benchtop onto paper toweling. Do not allow the wells to completely dry at any time.
  - 10.7.3 Add 300 µL of 1X Wash Buffer to each assay well.
  - 10.7.4 Incubate for 1 minute.
  - 10.7.5 Repeat steps 10.7.1 through 10.7.4 **two** more times.
- **10.8** Add 100 μL of prepared **1X Avidin-Biotin-Peroxidase Complex (ABC) into** each well, cover with plate sealer and incubate at 37°C for 30 minutes.
- **10.9** Wash plate 5 times with **1X Wash Buffer as** follows:
  - 10.9.1 Discard the liquid in the wells by rigorously flicking into an acceptable waste receptacle.
  - 10.9.2 Gently blot any remaining liquid from the wells by tapping inverted on the benchtop onto paper toweling. Do not allow the wells to completely dry at any time.
  - 10.9.3 Add 300 µL of 1X Wash Buffer to each assay well.
  - 10.9.4 Incubate for 1 minute.
  - 10.9.5 Repeat steps 10.9.1 through 10.9.4 four more times.
- **10.10** Add 90 μL of **TMB Color Developing Agent** to each well, cover with plate sealer and incubate at 37°C in the dark for 15-25 minutes.
  - (NOTE: optimal incubation time must be determined by the user. Optimal development can be visualized by blue shading in the top four standard wells, while the remaining standards are still clear.)
- **10.11** Add 100 μL of **TMB Stop Solution** to each well. Well color should change to yellow immediately.
- **10.12** Read the O.D. absorbance at 450 nm with a standard microplate reader within 30 minutes of stopping the reaction in step 10.11.



## 11. Calculation of Results

For analysis of the assay results, calculate the **Relative OD**<sub>450</sub> for each test or standard well as follows:

(Relative 
$$OD_{450}$$
) = (Well  $OD_{450}$ )–(Mean Blank Well  $OD_{450}$ )

The standard curve is generated by plotting the mean replicate **Relative OD**<sub>450</sub> of each standard serial dilution point vs. the respective standard concentration. The Human IGFBP7 concentration contained in the samples can be interpolated by using linear regression of each mean sample **Relative OD**<sub>450</sub> against the standard curve. This is best achieved using curve fitting software.


**Note:** if the samples measured were diluted, multiply the derived mean sample concentration by the dilution factor for a final sample concentration.

## 12. Typical Expected Data

**12.1** <u>Typical absorbance values.</u> Expected absorbance for standards when TMB incubation is performed for 20 minutes at 37°C and measured at OD<sub>450</sub>.

Standard Number	8	7	6	5	4	3	2	1
Standard Concentration (ng/mL)	0	0.625	1.25	2.5	5.0	10	20	40
OD <sub>450</sub>	0.042	0.147	0.213	0.413	0.809	1.227	1.849	2.347

**12.2** <u>Typical standard curve.</u> This standard curve is for demonstration purposes only. An assay specific standard curve should be performed with each assay.





# 12.3 Reproducibility

	Intra-Assay			Inter-Assay		
Sample ID	1	2	3	1	2	3
n	16	16	16	24	24	24
Mean Measured Concentration (pg/mL)	1629	4265	18882	1617	4427	20605
Standard Deviation	114.03	319.87	1340.62	118.04	358.58	1669
Consistency (%CV)	7	7.5	7.1	7.3	8.1	8.1



### 13. Technical Resources

#### 13.1 References

- 13.1.1 Murphy M, Pykett MJ, Harnish P, Zang KD, George DL (Dec 1993). "Identification and characterization of genes differentially expressed in meningiomas". Cell Growth Differ 4 (9): 715–22. PMID 7694637.
- 13.1.2 Yamauchi T, Umeda F, Masakado M, Isaji M, Mizushima S, Nawata H (Nov 1994). "Purification and molecular cloning of prostacyclin-stimulating factor from serum-free conditioned medium of human diploid fibroblast cells". Biochem J 303 (2): 591–8. PMC 1137368. PMID 7980422.
- 13.1.3 "Entrez Gene: IGFBP7 insulin-like growth factor binding protein 7".
- 13.1.4 Oh Y, Nagalla SR, Yamanaka Y, Kim HS, Wilson E, Rosenfeld RG (November 1996). "Synthesis and characterization of insulin-like growth factor-binding protein (IGFBP)-7. Recombinant human mac25 protein specifically binds IGF-I and -II". J. Biol. Chem. 271 (48): 30322–5.
- 13.1.5 Swisshelm, K., Ryan, K., Tsuchiya, K., Sager, R. Enhanced expression of an insulin growth factor-like binding protein (mac25) in senescent human mammary epithelial cells and induced expression with retinoic acid. Proc. Nat. Acad. Sci. 92: 4472-4476, 1995.

#### 13.2 Technical Support

#### **USA**

Aviva Systems Biology, Corp. 7700 Ronson Rd, Suite 100 San Diego, CA 92111

Phone: 858-552-6979 Toll Free: 888-880-0001 Fax: 858-552-6975

Technical support: techsupport@avivasysbio.com

#### China

Beijing AVIVA Systems Biology 6th Floor, B Building, Kaichi Tower #A-2 Jinfu Road. Daxing Industrial Development Zone Beijing, 102600, CHINA

Phone: (86)10-60214720 Fax: (86)10-60214722

E-mail: support@avivasysbio.com.cn

中国地址:北京大兴工业开发区金辅路甲 2 号凯驰大厦 B 座 6 层 (102600)

电话: 010-60214720/21 传真: 010-60214722

产品售前咨询及销售: sales@avivasysbio.com.cn

售后及技术支持: support@avivasysbio.com.cn